

ORTHODONTIC INSTRUMENTS



Table of Contents

Pin and Ligature Cutters	04	
Distal End Cutter		
Cinch Back	05	
Utility Pliers Weingart Pliers	06	
Bracket Remover Pliers	07	
Aligner Pliers	08	
Wire Bending Wire Forming Pliers	09	
Step Pliers	10	
Wire Bending	11	
	40	
Mathieu Needle Holders	12	
Gauge	13	
	4.4	
Hand Instruments	14	



















Pin and Ligature Cutters

- Wear and corrosion resistant steel jaws
- Steel jaws ensure consistent performance and reliability of cutting edges.
- Plier's shape with round edges allows maximum comfort and safety
- Forging is made with highest quality surgical grade stainless steel.
- Made with box-joint method, which stops lateral play or staggering.
- Retains cutting edges alignment and provides smooth and precise floating movements without effort.
- Exactly specified sizes and design of handles, ensures comfortable fit with the natural movement of hands.
- Stainless steel Matte finish with Mirror touch on the handle
- Narrow and slim jaws designed for easy access for hard to reach areas.
 Cuts soft wires, pins and elastics.





PR-1002

Micro Pin and Ligature Cutter Straight Soft wires up to 0.015"



PR-1006

Pin and Ligature Cutter Straight Soft wires up to 0.015"



PR-1004

Micro Pin and Ligature Cutter 15 Degree Angled Soft wires up to 0.015"



PR-1008

Micro Pin and Ligature Cutter 15 Degree Angled Soft wires up to 0.015"



PR-1018

Hard Wire Cutter Straight Archwire up To .022" X .028"



PR-1020

Hard Wire Cutter 15 Degree Angled Archwire up To .022" X .028"



Distal End Cutter/Flush Cutter





PROJSHARP ORWAL



PR-1022
Distal End Cutter Mini Handle
With Safety Hold
Maximum Cutting Capacity: .021" X .025"

PR-1026
Distal End Cutter Long Handle
With Safety Hold
Maximum Cutting Capacity: .021" X .025"

PR-1024
Distal End Cutter Standard Handle
With Safety Hold
Maximum Cutting Capacity: .021" X .025"







PRO/SHARP DENTAL

PR-1028
Flush Cut Distal End Cutter With Safety Hold Standard Handle
Maximum Cutting Capacity: .021" X .025"

PR-1030 Flush Cut Distal End Cutter With Safety Hold Long Handle Maximum Cutting Capacity: .021" X .025"

Cinch Back





Utility Pliers Weingart Pliers





PR-1032 Weingart Pliers Slim



PR-1036Weingart Pliers Angled Long Beak



PR-1038 Weingart Pliers Angled Short Beak



PR-1034 Weingart Pliers Utility



PR-1040 Weingart Pliers Large Large Jaw









Bracket Remover Pliers









PR-1050

Posterior Band Remover Pliers (Standard) Replaceable pads



PR-1060

Direct Bond Remover Pliers (Straight)

PR-1052

Posterior Band Remover Pliers (Long) Replaceable pads



PR-1062

Direct Bond Remover Pliers (Angled)

PR-1058

Adhesive Remover Pliers Replacement Carbide inserted Tip. Replaceable pads



Hole Punch

Aligner Pliers



Vertical Hole

Horizontal Hole



PRO/SHARP Wire Bending/Wire Forming Pliers

Standard Handle



PR-1072 Bird Beak Pliers Bends wire up to 0.030"





Bends wire up to 0.040"

PR-1074
Secure Hold Bird Beak Pliers Serrated Jaw



PR-1076
Long Beak Light Wire Pliers
Bends wire up to 0.020"



PR-1078 Micro Light Wire Pliers Bends wire up to 0.020"



PR-1082 Light Wire Pliers with Cutter 3-Grooved Bends & cuts wire up to 0.020"



PR-1090 Jarabak Pliers with Cutter Standard handle Bends & cuts wire up to 0.020"



PR-1092 Young Pliers Bends wire up to 0.040"



PR-1096 Adam Pliers Bends wire up to .040"




PR-1098 Lingual Arch Forming Pliers Bends wire up to 0.036"





PR-1120
De La Rosa Pliers Grooved
Grooved Size .016"- .018" - .022"



PR-1122 Ball Hook Crimping Pliers





PR-1124 Tweed Loop Forming Pliers Beak sizes .045" - .060" - .070" Wire up to 0.022"






Step Pliers







Hollow Chop Pliers Bends wire up to 0.036"

Optical Pliers

Nance Loop Pliers 3-4-5-6mm Bends wire up to 0.022"



Wire Bending



PR-1156 Ligature Wire Forming Pliers



PR-1140 Crown and Band Pliers




PR-1142 Band Contouring Pliers



PR-1144 Crown and Contouring Pliers



Wire up to 0.030"





PR-1094 Made with Stainless Steel. Easy to assemble



Mathieu Needle Holders









High Heat Bite Stick/Band Seater Autoclavable



All Purpose Scissor 4"





PRO/SHARP DENTAL

PR-1322 Direct Bond Bracket Holder With Fine Tip size: 5"



PR-1230 Mershon Band Pusher



PR-1318 Direct Bond Bracket Holder with Alignment Tip Size: 5"



PR-1322 Bracket Placing Tweezers with Alignment Tip Size: 5"

Gauge



PRO/SHAR

PR-1328 Molar Tube Holder Tweezers Reverse Action 5"



Boley Gauge PR-1232 Boone Gauge Stainless Steel Incisal Edge: 3.5, 4.0, 4.5, and 5.0mm



PR-1224 Positive Positioner





PR-1244 Bracket Height Gauge Movable Head



Hand Instruments



PR-1246 Ligature Director Curved Explorer Sharp



PR-1248 Ligature Director Straight/Angled



PR-1250 Medium Scalar Director



PR-1254 ligature taker double ended



PR-1298 Sickle Probe /Sharp Explorer



PR-1280 Spatula 6mm Single Sided



PR-1258 Band Pusher Small Scalar



PR-1264 Ligature Placing Double Ended



PR-1266Distal Bender
Double Ended



PR-1268Sickle Explorer
Angled Explorer



PR-1270 Sickle Explorer Single Sided



Angled Explorer Single Sided





Mouth Mirror with Handle PR-1286

Size: Size: #4 22mm

Instrument Sterilisation:

Inspect, clean, and sterilise before each use. Autoclave for at least 3 minutes at 134°C.

- Sterilisation is recommended before each use.
- Instruments must be cleaned and sterilised in the open position.
- Clean instruments thoroughly, rinse and dry before Sterilisation.
- Remove instruments immediately when cycle is complete to avoid corrosion or discoloration.
- Brown or yellow discolouration of instruments may occur if instruments are not dried after sterilisation, not well cleaned or not kept in a dry place.
- Discolouration may appear that may not be corrosion and can be removed by scrubbing with light wool saturate.
- Never Mix a corroded instrument with new instruments as oxidation may spread to new instruments.
- To help prevent corrosion, please use an alkaline ultrasonic detergent such as NCA 360.

PASSIVATION










HRC TESTING







PROSHARP DENTAL Unit 1, 4 Forge Place Narellan, NSW, 2567

Australia

www.prosharpdental.com.au

Email: Sales@prosharpdental.com.au

Phone: 1300 542 546